FR

Exacto®

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REF 0081N

1

BIOSYNEX
EASY DIAGNOSTICS FOR LIFE

SOMMAIRE

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AVITA Corporation
 9F, NO.78, Sec.1, Kwang-Fu Rd., san-Chung District 24158,
 New Taipei City,
 Taiwan, R.O.C.

SERVICE CONSOMMATEUR
 CONSUMER DEPARTMENT
 +33 (0) 3 88 78 85 24

client@biosynex.com
www.biosynex.com

Wellkang Ltd
 The Black Church, St. Mary's Place
 Dublin 7 , D07 P4AX , Ireland

20 Date de dernière modification : 10/2019

EN

Exacto® Blood Pressure wrist

REF 0081N

BIOSYNEX
EASY DIAGNOSTICS FOR LIFE

SUMMARY

Introduction	23
Important information before use	24 - 25
Blood Pressure Monitor Features	25
LCD Screen Description	26
Inserting Batteries	27
Setting the Date and Time	28
How to Put on Your Cuff	29
Measuring Heart Rate and Blood Pressure	30
Cardiac Arrhythmia Detector	31
World Health Organization (WHO)	31
Memory Function	32
Troubleshooting	33
Error codes	34
Care and maintenance	35
Standards Applied	35
Technical Specifications	36
EMC Tables	37 - 39

INTRODUCTION

Introduction

To get the maximum benefit from your blood pressure monitor, we recommend that you consult your doctor or qualified health care professional first.

Thank you for buying a blood pressure monitor. When properly used and carefully maintained, your blood pressure monitor will provide you with years of reliable readings.

The measurement method used by your automatic blood pressure monitor is called the oscillometric method. The blood pressure monitor detects your blood's movement in the artery in your wrist and converts the movements into a digital reading. The oscillometric method does not require a stethoscope, making the blood pressure monitor easier to use.

The pre-shaped cuff is easy to put on and comfortable to wear.

Clinical research has proven a direct relationship between blood pressure in the wrist and blood pressure in the arm. Because the arteries in your wrist and arm are

connected, changes in blood pressure in the wrist reflect the changes of blood pressure in the arm. Frequently measuring blood pressure in the wrist will provide you and your doctor with an accurate indication of changes in your actual blood pressure.

Those with severe hypertension, severe arteriosclerosis or severe diabetes must consult their doctor regarding blood pressure monitoring in the arm.

Blood pressure measurements with this device are equivalent to measurements obtained by a trained observer using the cuff/stethoscope auscultation method within the limits prescribed by the National Standard for Electronic or Automated Sphygmomanometers.

IMPORTANT INFORMATION BEFORE USE

When you use this product, be sure to follow all instructions below. Any action contrary to these instructions may affect the accuracy of the device.

1. Do not modify this equipment without the authorization of the manufacturer.

2. Blood pressure measurements should only be interpreted by a qualified medical professional who knows your medical history.

By using this device regularly and recording your measurements, you can keep your doctor informed of changes in your blood pressure

3. Perform your measurements in a quiet place.

You should be seated in a relaxed position.

4. Avoid smoking, eating, taking medication, consuming alcohol, and any physical activity 30 minutes before taking your blood pressure.

If you are showing signs of stress, avoid taking your

blood pressure until the feeling has faded.

5. Rest for 15 minutes before taking your blood pressure.

6. Take off tight clothing or jewellery that may interfere with the placement of the cuff.

7. Keep the blood pressure monitor stable while measuring to get an accurate reading.

Keep still and do not speak during the measurement.

8. Record your daily blood pressure and pulse readings.


9. Always take your blood pressure at the same time each day or according to your doctor's advice to get an accurate indication of your actual changes in blood pressure.

IMPORTANT INFORMATION BEFORE USE


- 10.** Wait at least 15 minutes between measurements to allow blood vessels to return to normal.
- 11.** In the event that your pulse is very weak or irregular, errors may occur that prevent correct measurement. If any abnormal changes are observed, consult your doctor or health care professional.
- 12.** This device is intended for adult use. When taking a measurement, you can stop the cuff inflation or deflation at any time by pressing the POWER button.
- 13.** If you need help during the setup, use or care, please contact the manufacturer.

Using a blood pressure monitor cannot in any case exempt you from consulting a doctor. Do not change your treatment without medical advice.

BLOOD PRESSURE MONITOR FEATURES

LCD SCREEN DESCRIPTION


INSERTING BATTERIES

Low battery warning:

You must replace the batteries when the low battery icon «  » appears on the screen or when the screen does not turn on after pressing the POWER button.

To change the batteries:

1. Open the battery compartment cover.
2. Make sure the polarity of the batteries is correct and insert two AAA batteries
3. Close the battery compartment cover.


Recommendations

1. Always use new batteries; never mix new and old batteries.
2. Do not use rechargeable batteries; rechargeable batteries may vary in quality and power. The use of rechargeable batteries may affect the performance of this device.
3. Removing the batteries is recommended if the device is not being used for a prolonged period of time.
4. Batteries must be disposed of in accordance with environmental standards and institutional policies.

SETTING THE DATE AND TIME

You must set the device's date and time before the first use and each time the batteries are replaced

1. When the device is turned off (blank screen) or after replacing the batteries, press and release the «Date and time» button; «year» will begin to flash on the screen.
2. Press the «MEMORY» button to display the desired year, press on «POWER» to confirm the «Year».
3. The «Month» will flash next. Repeat steps 1 and 2 to set the «Month» and the «Date», then «Hour», then «Minute».
4. After setting the minute, the device will automatically exit the date/time setting and quickly show the word OFF before turning off.


HOW TO PUT ON YOUR CUFF

It is important to avoid smoking, eating, taking medication, consuming alcohol, and any physical activity 30 minutes before taking your blood pressure. If for any reason, you cannot or should not use your left wrist, modify the instructions to apply the cuff on your right wrist. Your doctor can help you identify the best wrist for taking the measurement.

1. Remove any tight clothing or jewellery that may interfere with positioning the cuff.
2. Sit at a table or desk with your feet flat on the floor.
3. Place the back of your hand on the blood pressure monitor box, palm up, so that the device is at the same level as your heart.
4. Apply the pre-shaped cuff to your wrist so the digital display is facing you.
5. Adjust the cuff to about 1 cm away from your hand.
6. The cuff should fit comfortably, but be snug around your wrist.

Blood pressure varies naturally from one wrist to another; measure your blood pressure on the same wrist so that the measurements are comparable.

MEASURING HEART RATE AND BLOOD PRESSURE

Please read this entire manual before taking your first measurement.

1. Rest your elbow on a hard surface with your palm up. Raise your wrist so the blood pressure monitor is at the same level as your heart. Relax the hand of your other arm.

2. Press the POWER button to start the device.

3. After the automatic verification, the values of the last measurement will appear on the screen.

4. The blood pressure monitor will begin the measurement.

5. The cuff will inflate automatically, the screen indicating the increasing pressure in the cuff. When the pressure increases, an arrow pointing up will appear on the screen.

6. When the inflation has reached the optimal level, the screen will display the decreasing pressure, showing an arrow pointing down, and you will feel the pressure in the cuff decrease.

7. To detect the heart rate, the heart rate icon will appear and flash continually on the LCD screen.

8. Your blood pressure and pulse measurements will show simultaneously on the screen.

9. The hypertension indicator will indicate your reading range separately on the screen.

10. The measurement is automatically stored in the memory.

11. Select POWER to turn the device off and save energy, battery life. The device will turn off automatically after two minutes.


IRREGULAR HEARTBEAT DETECTOR

1. Your digital blood pressure monitor has an irregular heartbeat detector. An irregular heartbeat may influence the measurement. If the blood pressure monitor detects an irregular heartbeat during the measurement, the «  » icon will appear on the screen with the measurement values. You may take another measurement to ensure that the values are not being influenced by movement during the reading or by an irregular heartbeat.


2. If this icon appears frequently, we advise you to consult with your doctor.

3. This device does not replace a cardiac exam, but serves to detect heartbeat irregularities at an early stage.

WORLD HEALTH ORGANIZATION (WHO)

This device is equipped with our unique hypertension indicator. The World Health Organization has established globally recognized standards to evaluate high or low blood pressure readings. The graphic below should only be considered a guideline; always consult your doctor or health care professional to interpret your individual results.


MEMORY FUNCTION

Recall measurements in the Memory:

You can recall up to 60 total measurements in the memory, plus an average of the 3 last measurements saved in the memory to share them with your doctor or a qualified health care professional.

1. Press and release the «MEMORY» button. The device will then show the average of the last 3 saved measurements.
2. Continue to push the «MEMORY» button to see the previously saved measurements. The measurements will appear on the screen from the most recent to the oldest; the number of memory will appear on the lower left corner.

Press the
«MEMORY»
button.

3. All the results of a given measurement will show, including the measurement results, the heart rate, the hypertension indicator and the timestamp.

4. When the number of readings surpasses 60, the oldest data will be replaced by the new recording.
5. Press the Power button to turn off the blood pressure monitor at any time while reading saved measurements

Erasing Measurements from the Memory:

From the turned off screen, hold the «MEMORY» button down until the screen shows CLr. This indicates that all measurements have been erased.

TROUBLESHOOTING

Issue	Likely cause	Recommended action
Nothing appears on the screen even when the device is powered on.	The batteries are dead.	Replace all batteries with new batteries.
	The batteries are not properly aligned with the tracks.	Replace the batteries in the proper position.
The low battery icon appears.	The batteries are dead.	Replace all batteries with new batteries.
	In very cold conditions, the batteries have lower electrical charges.	Warm up the batteries or use the device in a warmer environment.
The device's operating time is inconsistent.	Different brands of batteries have different lifetimes.	Use alkaline batteries and replace all the batteries at the same time with the same brand of batteries.
No reading after the measurement.	The batteries are dead.	Replace all batteries with new batteries.
Questionable blood pressure results.	The cuff may have been poorly positioned.	Adjust the patient and the wrist cuff to measure.
	Blood pressure varies naturally throughout the day.	Rest for a moment, relax and measure again.
Questionable heart rate results.	Body movement during use of the device.	Avoid moving while taking the reading.
	Measurements taken shortly after exercise or time outside.	Do not take measurements after exercise or coming back in from outside.
The device turns off automatically.	Device design.	Press the power button again, then measure again.
During the measurement, the cuff begins to inflate again.	This could be a normal action if the user's blood pressure is higher than the initial pressure value, the device automatically pumps to a higher pressure level of 40mmHg each time.	Relax, and try to take a measurement again.
	The wrist cuff is not secured.	Make sure that the wrist cuff is securely fastened and retake the measurement.

ERROR CODES

Error code	Meaning	Corrective action
Err 0	No pulse or insufficient pulse detection	Remove heavy clothing and retry.
Err 1	Leak in the pressure cuff/inflation too weak	The wrist cuff is not securely fastened. Put it on again and retake the measurement.
Err 2	Pressure failure	Rest for a moment, relax and retry.
Err 3	Deflation failure	The wrist cuff is not securely fastened. Put it on again and retake the measurement.
Err	Memory error	Take out the batteries to restart the device, then retake the measurement.
	Low batteries	Replace all batteries with new batteries.

CARE AND MAINTENANCE

- 1.** Clean the device and wrist cuff carefully with a slightly damp cloth.
- 2.** Do not put the device in water. Damage could occur if water gets into the device.
- 3.** Do not use gas, powerful detergent, or solvents to clean the device, including the wrist cuff.
- 4.** Use a soft dry cloth to clean your blood pressure monitor.
- 5.** Do not store in direct sunlight, dust, or excessive humidity.
- 6.** Avoid extreme temperatures.
- 7.** Do not dismantle the blood pressure monitor or the cuff.
- 8.** Remove the batteries if the blood pressure monitor is not being used for a prolonged period of time.

STANDARDS APPLIED

This product is in accordance with the measures in directive CE MDD (93/42/CEE). The following standards apply to the design and/or the manufacturing of the products:

EN 1060-3

Non-invasive sphygmomanometers – Part 3: Additional requirements for the electromechanical blood pressure measuring system

EN 1060-4

Non-invasive sphygmomanometers – Part 4: Test procedures that aim to determine the overall system accuracy of automated non-invasive sphygmomanometers

ISO 14971

Medical Devices – Application of risk management to medical devices. Classification according to IEC/EN 60601-1 paragraph 5:

- Internally-powered equipment
- IPX0
- Equipment not suitable for use in the presence of an anaesthetic mixture flammable with air or oxygen or nitrous oxide.
- Continuous operation

TECHNICAL SPECIFICATIONS

• Measuring range :	Blood pressure: 30 to 280 mmHg Pulse rate: 40 to 199 beats/min
• Calibration accuracy :	Blood pressure: \pm 3 mmHg Pulse rate: \pm 4% of reading
• Operational environment:	10°C to 40°C With relative humidity up to 85% (without condensation)
• Storage/transportation environment:	-20 ° C to + 50 ° C With relative humidity up to 85% (without condensation)
• Power :	2 AAA 1,5 V Batteries (LR03)
• Weight :	approx. 120g (excluding batteries)
• Dimensions :	approx. 71 mm x 84 mm x 32.2 mm (L x H x D)
• Cuff circumference (size M):	approx. 12.5 ~ 21 cm

EMC TABLES

Manufacturer's declaration and directive – electromagnetic emissions

The 0081N wrist blood pressure monitor is intended for use in the electromagnetic environment specified below. The client or user of the 0081N wrist blood pressure monitor must be sure to use it in such an environment.

Emissions test	Compliance	Electromagnetic environment guidelines
RF Emissions CISPR 11	Group 1	The 0081N wrist blood pressure monitor uses RF energy only for its internal operation. Therefore, its RF emissions are very low and are not likely to cause interference in electronic equipment nearby.
RF Emissions CISPR 11	Class B	
Harmonic Emissions CEI 61000-3-2	Not applicable	
Voltage variations / Flicker emissions CEI 61000-3-3	Not applicable	The 0081N wrist blood pressure monitor may be used in all establishments, including dwellings and establishments directly connected to the public low-voltage power supply system that powers buildings intended for homes.

EMC TABLES

Manufacturer's declaration and directive – electromagnetic immunity

The 0081N wrist blood pressure monitor is intended for use in the electromagnetic environment specified below. The client or user of the 0081N wrist blood pressure monitor must be sure to use it in such an environment.

Immunity Test	CEI 60601 Test level	Compliance level	Guidelines on the electromagnetic environment
Electrostatic Discharge (DES) CEI 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	The floors must be made of wood, cement or ceramic tiles. If the floors are covered in synthetic material, relative humidity must be below 30%.
Power Supply Frequency (50/60 Hz) Magnetic Field CEI 61000-4-8	3 A/m	3 A/m	The magnetic field of the supply frequency must be at a level appropriate to the place of use.

EMC TABLES

Manufacturer's declaration and directive – electromagnetic immunity

The 0081N wrist blood pressure monitor is intended for use in the electromagnetic environment specified below. The client or user of the 0081N wrist blood pressure monitor must be sure to use it in such an environment.

Immunity Test	CEI 60601 Test level	Compliance level	Electromagnetic environment guidelines
Radiated RF CEI 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	<p>Portable and mobile RF communication equipment must not be used closer to any part of the 0081N wrist blood pressure monitor, including the cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance $D = 1.2 \sqrt{P}$ $D = 1.2 \sqrt{P} \text{ 80 MHz to 800 MHz}$ $D = 2.3 \sqrt{P} \text{ 800 MHz to 2.5 GHz}$ Where P is the maximum output power of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>The field strengths of fixed RF transmitters, as determined by an electromagnetic site survey, must be less than the compliance level in each frequency range.</p> <p>Interference may occur near equipment marked with the following symbol: </p>

NOTE 1 At 80 MHz and 800 MHz, the highest frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by the absorption and reflection of structures, objects and people.

A. The field strengths of stationary transmitters, such as base stations for radio (cellular/wireless) and land mobile radios, amateur radios, AM and FM broadcasting and TV broadcasting may not in theory be predicted accurately. To evaluate the electromagnetic environment due to fixed RF transmitters, consideration should be given to an electromagnetic survey of the site. If the measured field strength value at the location of the 0081N wrist blood pressure monitor exceeds the applicable RF compliance level above, the 0081N wrist blood pressure monitor should be checked to verify normal operation. If abnormal performance is observed, additional measures may be required, such as reorienting or reinstalling the 0081N wrist blood pressure monitor.

B. Depending on the frequency range from 150 kHz to 80 MHz, the field strengths must be less than 3 V/m.

CEM TABLES

Recommended separation distance between portable and mobile RF communication equipment and the 0081N wrist blood pressure monitor.

The 0081N wrist blood pressure monitor is intended for use in an electromagnetic environment where the RF perturbations are controlled. The customer or user of the 0081N wrist blood pressure monitor can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications (Transmitters) and the 0081N wrist blood pressure monitor as recommended below, depending on the maximum output power of the communication equipment.

Maximum rated output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2,5 GHz
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters whose maximum output power is not specified above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the Maximum output power of the transmitter in watts (W) according to the manufacturer of the transmitter.

NOTE 1 At 80 MHz and 800 MHz, separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by the absorption and reflection of structures, objects and people.

	Read the instruction manual before use
	Type BF applied part
	Batteries and electronic devices must be disposed of in accordance with the locally applicable regulations
	Manufacturer
	Storage conditions
	Catalog number
	Keep dry
	Lot number
	This device complies with the requirements of the Medical Device Directive 93/42/EEC.
	Production year
	Common signage TRIMAN
	European Authorized Representative

WARRANTY

This product is warranted for a two-year period from the purchase date. However, this warranty does not apply to damages, in particular, resulting from non-compliance with this guide, incorrect manipulations, shocks or inappropriate installation of the battery. It does not apply either to the battery or the package.

Distribué par:

BIOSYNEX SA
22 Boulevard Sébastien Brant
67400 Illkirch-Graffenstaden
FRANCE

AViTA Corporation

9F, NO.78, Sec.1, Kwang-Fu Rd., san-Chung District 24158,
New Taipei City,
Taiwan, R.O.C.

SERVICE CONSUMERS
CONSUMER DEPARTMENT
+33 (0) 3 86 78 85 24

client@biosynex.com

www.biosynex.com

Wellkang Ltd
The Black Church, St. Mary's Place
Dublin 7 , D07 P4AX , Ireland

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AU S.A.V. BIOSYNEX L'ADRESSE INDIQUEE CI-DESSOUS