# TEST TETANOS

Test destiné à l'auto-diagnostic

# Exacto®

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


# TETANUS TEST

For self-testing

**EXACTO® Tetanus Test** is an easy-to-use, reliable self-test kit that can be used at home to determine tetanus vaccine status.

Tetanus is a serious infection caused by the neurotoxin produced by the Clostridium tetani bacterium. This bacterium, naturally present in soil and animal excrement, is highly resistant and can enter the body through a simple wound. Once in the body, the bacterium releases the tetanus toxin which blocks the release of neurotransmitters by neurons and causes muscle spasms and muscle paralysis.

Vaccination remains the only way of preventing *C. tetani* infections in children and adults. Immunoprotection against the tetanus toxin is assessed by measuring the levels of anti-tetanus toxin antibodies in the blood.

Knowing your tetanus antibody level helps to determine if you are effectively protected against Clostridium tetani infections. A negative result with the **EXACTO® Tetanus Test** self-test indicates an insufficient level of protection.

**EXACTO® Tetanus Test** is a rapid immunochromatographic test capable of detecting anti-tetanus toxin antibodies in one drop of blood, enabling you to find out your tetanus vaccine status. In the event of a negative result, a vaccination booster will be sufficient to ensure protection against tetanus.

## KIT CONTENTS

## TIPS FOR DRAWING BLOOD

Wash your hands in warm water beforehand and dry them properly. Warm water improves circulation in the fingertips.

You can let your hand hang down by your side for about a minute so that more blood flows into it. We recommend that right-handed people prick their left hand and left-handed people prick their right hand.

Preferably prick the side of your fingertip: this is where there is the most blood and the least sensitivity. Avoid areas where the skin may be thicker and more difficult to puncture.

It may be that a drop does not form immediately. In this case, you should wait a little while and massage your finger gently, from the base to the puncture site.

Feel free to form large, hanging blood drops at the fingertip. This will make it easier and quicker to fill the pipette.

## DRAWING

**1** Please read the directions for use in their entirety. Make sure you have a watch or stopwatch to hand (the only equipment not provided).

**2** Wash your hands with soap and warm water, then dry them before moving on to the next step.

**3** Open the box and **remove all components from their packaging**. Place them on a clean, flat surface. Identify each item of the kit (see kit contents on the left, and make sure you have all 7 items, A to H).

Tear along the notch in the pouch **A**. Take out the cassette **B** and dispose of the desiccant. **Use the test within one hour.**

Unscrew the cap on the bottle of diluent **C** and place it upright on the work surface.


Take the alcohol wipe **E** out of its case and disinfect your fingertip. We recommend that you prick the middle finger (see photo).


Take the compress **F** out of its case and wipe off any traces of alcohol.


Pull the cap off the lancing device **G**.

Place the lancing device **G** against the fingertip and press firmly to release the lancet.

**Squeeze the pricked fingertip firmly to form a large drop of hanging blood.**

Place the pipette **D** in direct contact with the drop of blood.

Fill the pipette **D** with blood completely, up to the white filter.

## TESTING

Empty the contents of the pipette into the round well **S** in the cassette **B** by **pressing down the white piston of the pipette.**

Place **3 drops of the diluent C** into the round well **S** of the cassette **B**.

Start the stopwatch or note the time this is done. **Wait for 10 minutes before reading the result.** The test should not be read after 15 minutes. In the meantime, you can put on the dressing **H**.

## INTERPRETATION

**NEGATIVE**

**POSITIVE**

**INVALID**

**PUT YOUR TEST HERE**

The test result is **NEGATIVE** if a coloured line appears only in the control zone (C). No coloured line appears in the test zone (T). A negative result indicates the absence of anti-tetanus toxin antibodies and therefore no protection against tetanus. **Consult your doctor as soon as possible!**

The test result is **POSITIVE** if two coloured lines appear:  
- one line in the control zone (C)  
- one line in the test zone (T).  
A positive result indicates the presence of sufficient quantities of anti-tetanus toxin antibodies. **You are effectively protected against tetanus.**

The result of the test is **INVALID** if no coloured line is visible in either zone, or if no coloured line appears in the control zone (C). **You should perform another test.**

## PERFORMANCE

The **EXACTO® Tetanus Test** has formed the subject of various studies. A 91% correlation was observed between the **EXACTO® Tetanus Test** and the ELISA assay.

- A clinical study conducted on 782 serum samples showed that the **EXACTO® Tetanus Test** has a sensitivity of 96% and a specificity of 88%.
- A practicability study conducted in a lay population showed that 100% of participants obtained a correct and interpretable result.
- A study on the reading and interpretation of the test confirmed that 100% of the various types of results were interpreted correctly.

## FAQ

**What is the purpose of the control line (C) ?**  
The appearance of a line in the control zone (C) of the cassette is an internal migration check, indicating that the test has functioned correctly.

**A line appears in the test zone (T) of my EXACTO® Tetanus Test but it is faint. Is my result really positive ?**  
The intensity of the line does not affect the reading of the results. As such, any line that appears in the test zone (T) must be interpreted as a positive result, regardless of its intensity.

**What do I do if the result of my EXACTO® Tetanus Test is negative ?**  
A negative result with the EXACTO® Tetanus Test means that you are not protected against the tetanus toxin. Consult your general practitioner as soon as possible to review your vaccination schedule. If the test is performed immediately following an injury with a risk of infection and the test result is negative, go to the nearest Emergency Department to receive appropriate care.

**What do I do if the result of my EXACTO® Tetanus Test is positive ?**  
A positive result with the EXACTO® Tetanus Test means that sufficient quantities of anti-tetanus toxin antibodies are present in your blood. You are therefore effectively protected against tetanus. Nevertheless, it is advisable to check when your vaccination boosters are due. Vaccination against tetanus is mandatory for infants and boosters are recommended in childhood, adolescence and adulthood. Ask your general practitioner for advice.

## PRECAUTIONS FOR USE

- This test must not be frozen and must be stored in a dry place at between 2°C and 30°C.
- Remove the test cassette from the airtight pouch and use it within one hour.
- This test must be used at a temperature of between 15°C and 30°C.
- This test must be used immediately after the drop of blood has been drawn from the fingertip.
- This test is for single use only.
- Make sure that the blood has been placed in the well before adding the drops of diluent.
- Keep the self-test kit and its components out of the reach of children; the buffer solution can present a danger in the case of ingestion and the lancet a risk of injury.
- Do not use the self-test after the expiry date given on the aluminium foil pouch and on the outside of the cardboard box.
- If the result is negative, consult your doctor to review your vaccination schedule.
- This test is not a substitute for a medical consultation or the results of a biological analysis carried out in a medical analysis laboratory.
- Erroneous results can be obtained if:
  - the test has not been used in accordance with the directions for use,
  - the foil pouch is damaged or if the test is not performed immediately after opening the foil pouch,
  - the storage conditions are not adhered to or if the test is performed after the expiry date indicated on the foil pouch.
- The used test cassette and all test components may be disposed of with household waste.
- All components included in this kit are intended for this test only. Do not reuse this test and/or its components.

Caution, see the directions for use	Number of tests per kit	Manufacturer	Lancet	REF 858963
in vitro diagnostic medical device	Use before	Do not re-use	TRIMAN unified signage	Pipette
Diluent	Batch Number	Catalogue No.	Green point	

Date of last revision : 10/2019